

POST-MARKETING STUDY

OBJECTIVES:

The goal of this prospective, open-label study was to demonstrate the safety and effectiveness of this device for reducing subcutaneous abdominal fat using an improved power delivery curve.

METHODS:

Male and female subjects with Fitzpatrick skin types I-VI ($N=26$) were treated. Up to four abdominal zones, up to 150cm² each, customized in size and location for body habitus were treated. Each zone underwent a single 20-minute treatment session. Follow-up visits occurred after 6 and 12 weeks. Using a standardized protocol, ultrasound measurement of subcutaneous abdominal fat thickness, abdominal circumference, reported patient satisfaction and digital images were obtained.

RESULTS:

The mean treatment area was 378.5 cm². At Week 12, there was a 21.6% mean reduction in abdominal fat thickness and a 1.6-inch mean reduction in abdominal circumference. Most subjects (96%) were satisfied or very satisfied with their results. The mean pain score was 2.5 on an 11-point ordinal scale. There were no non-responders. Two adverse events were mild transient erythema ($n=1$) and localized subcutaneous firmness ($n=1$) which resolved spontaneously within 12 weeks.

| Study | Dominion 2018 FDA Study | Dominion 2019 Post-Marketing Study |
|-----------------------------|-------------------------|------------------------------------|
| Sample Size | 36 | 26 |
| Treatment Time | 15 minutes | 20 minutes |
| Week 12 Mean Fat Reduction | -3.4 mm | -6.3 mm |
| Percentage of Fat Reduction | -15.1% | -25.3% Lower Abdomen |
| Subject Satisfaction | 89% | 96% |
| % Subjects with Nodules | 11.1% (4/36) | 3.8% (1/26) |

CONCLUSION:

EON is safe and effective for reducing abdominal fat and represents an improvement on the prior treatment protocol.